



Click-Tip™ Injection Needles

Designed for controlled needle deployment and injection, CONMED offers a full line of sclerotherapy needles to meet your needs.

CONMED's Click-Tip™ needle family includes 25, 22 and 19 gauge needles for accommodating a variety of clinical situations. The ergonomic, nurse-designed handle is easy to use for maximum comfort.



Features and Benefits

Top Down Injection

Allows for easy injection of solutions while the parallel port enhances handling and injection ergonomics.

Multiple Gauges and Lengths

CONMED has a full line of needles to meet your needs and accommodate a variety of procedure types.

Color Indication on "Unlock Button"

Clearly shows how to unlock needle in a low-light setting.

Needle Control System Designed for Safety

Safety Tab decreases the risk of scope damage and limits risk of sharps injury while the locking feature helps minimize needle bounce.

15° Beveled Stainless Steel Needle

Provides maximum piercing capabilities for easy insertion.



**To arrange a demonstration, call 800-448-6506.
To locate your local CONMED Sales Representative
or visit CONMED.com for more information.**

To learn more about this and other innovative products, call
1-800-448-6506 or visit **CONMED.com** for more information.

Ordering Information

Description	Needle Length	Catheter O.D.	Catheter Length	Qty	Catalog Number
CLICK-TIP™ Injection Needle					
Click-Tip™ Injection Needle, 19 Gauge	4mm	2.3mm	180cm	10/box	02-23-180
Click-Tip™ Injection Needle, 19 Gauge	4mm	2.3mm	230cm	10/box	02-23-230
Click-Tip™ Injection Needle, 19 Gauge	6mm	2.3mm	180cm	10/box	04-23-180
Click-Tip™ Injection Needle, 19 Gauge	6mm	2.3mm	230cm	10/box	04-23-230
Click-Tip™ Injection Needle, 22 Gauge	4mm	1.9mm	180cm	10/box	01-19-180
Click-Tip™ Injection Needle, 22 Gauge	4mm	1.9mm	230cm	10/box	01-19-230
Click-Tip™ Injection Needle, 22 Gauge	6mm	1.9mm	180cm	10/box	03-19-180
Click-Tip™ Injection Needle, 22 Gauge	6mm	1.9mm	230cm	10/box	03-19-230
Click-Tip™ Injection Needle, 25 Gauge	4mm	1.9mm	180cm	10/box	05-19-180
Click-Tip™ Injection Needle, 25 Gauge	4mm	1.9mm	230cm	10/box	05-19-230
Click-Tip™ Injection Needle, 25 Gauge	6mm	1.9mm	180cm	10/box	06-19-180
Click-Tip™ Injection Needle, 25 Gauge	6mm	1.9mm	230cm	10/box	06-19-230

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com